



Because others look to you...
Look to VisionTek®

Slide review with live digital microscopy

 **VisionTek®**
SAKURA
Digital Microscope
Brightfield System



GAIN A NEW PERSPECTIVE INNOVATIVE 3-IN-1 MICROSCOPY

Live diagnosis

VisionTek® is a desktop digital pathology microscope that offers live viewing with remote control for digital image capture of anatomical pathology and cytology slides. It is a fully motorized brightfield optical system with high-resolution fast cameras and a rapid autofocus to capture high-resolution live images.

- Review slides side-by-side, and navigate, magnify, measure, and annotate live
- Control your view remotely via the Internet, and review and share live images with multiple pathologists
- View live images ergonomically

Remote consultation


VisionTek®, coupled with broadband Internet access, enables live consultations anywhere.

- Frozen section slides can be reviewed remotely without a trip to the surgical suite
- Rural hospitals can have instant access to pathology specialists
- Clinicians can gain access to specialists worldwide for immediate expert diagnostic review

Scan archiving

VisionTek® is capable of user-defined Selective Slide Image (SSI™) scanning of multiple z-stacks, automated whole-slide scanning, or Partial Slide Image (PSI™) scanning by auto tissue detection.

- Diagnose live and then archive only those images you deem valuable
- Scan selective or partial scans for generating reports
- Scan high-resolution images for further analysis or educational purposes



Digital Microscope Brightfield System

Case review

VisionTek® enables complex cases to be reviewed **live** or in remote conference by multiple pathologists.

- *Share pathology findings with clinicians in tumor boards*
- *Compare and review all staining methods side-by-side*
- *Live or stored images can be shared for consultation*



Diagnosis

VisionTek® facilitates management of your workload and fine-tunes your review of complex cases.

- *The Virtual Case™ feature enables easy viewing, capture, and management of image data from multiple slides of individual cases*
- *Allows you to compare and review all staining methods side-by-side*
- *Facilitates live and remote examination and consultation of slides including frozen sections*



Research

VisionTek® can be your workhorse in a research lab where routine workflow mimics the clinical environment.

- *Allows rapid sorting of slides that need to be scanned for further analysis, archiving, and publication*
- *Enables rapid viewing and comparison of tissue microarrays (compare up to 16 images on the screen)*
- *Facilitates remote sharing of live images with your research collaborators*



Education

VisionTek® offers a versatile alternative to your standard multiheaded microscope.

- *Share images live on a large overhead monitor for viewing by an entire study group*
- *Enable remote study groups to control and review images when and how they choose*
- *Measure and annotate live images for a dynamic teaching experience*



VisionTek® DIGITAL MICROSCOPE/BRIGHTFIELD SYSTEM

THE VIEW HAS NEVER BEEN BETTER

Instrument Dimensions		Optics	
Height:	45.5 cm (17.9 in)	Magnification:	Slide tray overview 0.45x
Width:	40.5 cm (15.9 in)	Live image:	2.5x; 10x; 20x; 40x (digital zoom)
Depth:	52.0 cm (20.5 in)	Illumination:	White LED, with Kohler illumination
Weight:	35 kg (77.2 lb)	Z-stage autofocus:	Mechanical high-precision autofocus; motor controlled via software (mouse); IR beam and image-based autofocus mechanisms work together to focus automatically after each XY movement
		XY-stage automation:	Accommodates 4 standard microscope slides on a slide carrier; controlled via software
Image Detection		Computer	
Overview camera:	Resolution at 0.45x = 12 $\mu\text{m}/\text{pix}$	Operating system:	Windows 7, 64
Live view and scanning camera:	Resolution at 2.5x = 2.2 $\mu\text{m}/\text{pix}$ at 10x = 0.55 $\mu\text{m}/\text{pix}$ at 20x = 0.275 $\mu\text{m}/\text{pix}$	Processor:	Intel Core i7 or higher
Scanning time:	3 min for 15 mm x 15 mm (at 0.275 $\mu\text{m}/\text{pix}$)	Hard disk:	128 GB + 1000 GB 3.5"; 1000 Base-T Ethernet Intel, Pro/1000PT server network card
		Monitor:	24-inch 1920 x 1200 pixel resolution; display port or DVI

Operation Modes

Microscope Mode

- Overview of all slides
- Live review of single or multiple slides
- Scanned image review of single or multiple slides
- Single snapshot or user-defined slide scan
- Annotation and measurement

Scanner Mode

- Auto tissue detection mode
- Predefined scan mode

Digital Slide Gallery

- Manage and review scanned images

Visit our web site at www.sakuraus.com

The snapshot, partial scan, and whole slide scanned images are for research use only (RUO) in the USA and may not be used in diagnostic procedures.

© 2013 Sakura Finetek USA, Inc. All rights reserved.

0007000-01 Rev. A.



Sakura Finetek USA, Inc.
1750 West 214th Street
Torrance, CA 90501 USA
Phone: (800) 725-7823
www.sakuraus.com

Printed in USA